

Monitoring immunogenicity of your candidate vaccine by ELISpot

How to reduce assay variability?

Alain Poyau, Claire Serraz, Audrey Pabiou, Aurélie Chene, Lucie Levast, Marc Auvergnon, Luc Boban, Florian Brajon, Marie-Ambre Monet, Romain Roth, Caroline Seiner, Julie Silva, Cécile Valverde, Cendrine Josson, Christine Bain-Wendlinger

1. KCAS Bio www.kcasbio.com | +33 (0) 4 37 70 87 00 | 60F Rockefeller Av., Bioserra II – 69008 – Lyon – France

UNDERSTANDING THE CONTEXT OF USE OF THE ASSAY

Understanding Science



- What's the indication?
- What is the expected magnitude of responses?
- What's the target of your product?
- What is the expected impact of your product on the immune system?

Organizational constraints of the clinical study

- Starting date & duration of the study
- How many clinical sites & where are they located?
- What is the injection & sampling schedule?
- Who are the other stakeholders? Is there a central laboratory, a CRO, a logistical partner....?





and what will the data be used for?



- Monitor changes from baseline?
- Compare cohorts of patients?
- Is it a primary, secondary or exploratory endpoint?
- Make decisions on route, schedule, formulation, dose of injection...?
- Regulatory submission?

HIGH QUALITY SAMPLES

SHIFTING THE PARADIGM FROM CENTRALIZATION TO MULTI-SITE LOCAL LABORATORIES



Central Lab

- Standard SOP
- Samples from 28 sites in US & EU
- Samples shipped overnight

Clinical sites (n=3)

- Non-standardized protocol
- PBMC prepared and frozen within 8h from venipuncture

PBMC network (ACTIVE)

- Laboratory audited, trained and qualified
- Standardized protocol
- PBMC prepared and frozen within 8h from venipuncture

METHOD PERFORMANCE

I. IFN-γ ELISpot METHOD THAT FITS YOUR NEEDS

A. Antigen	B. Peptide conc.	C. Standardization	D. Assay format / target responding cells
Recombinant NP	Conc A <input checked="" type="checkbox"/>	Human plasma	Ex-vivo (circulating effector T cells)
Peptide pool <input checked="" type="checkbox"/>	Conc B	Serum-free medium <input checked="" type="checkbox"/>	After 5-12-days <i>In vitro</i> expansion (circulating central memory T cells)

: retained condition for clinical runs

PBMC

- From 3 healthy donors (A, D)
- From 6 patients (B)


Stimulation

- HCMV pp65 peptide pool (PP) (D); Test PPI-2 (A, D); Individual peptides from test PP (P1-3); CEFT peptide pool (C) Negative control (Med)
- Recombinant target antigen (Rec Ag)

IFN-γ ELISpot

- Ex vivo (A, D)
- 5-days expansion (B, C, D)
- 10-days expansion (D)

II. VARIABILITY OF THE METHOD DEPENDS ON ELISpot DESIGN AND LEVELS OF RESPONSES



Samples:

- In-study samples with different levels of response

Antigen:

- Vaccine-derived peptide pool

ELISpot:


- Ex-vivo
- After 5-days expansion

Design:

- 3 independent runs
- At 3 different days
- By 3 operators

III. CONTROL CHART & QUALITY MONITORING

A. Controlling suitability of the method along clinical testing



Samples:

- Reference sample (QC)
- Antigen: CEFT peptide pool

Design:

- Acceptance limits = mean +/- 3SD of first 10 analytical runs
- Over an 8-month period
- 3 operators

B. External Quality Assessment

- Proficiency testing with EQAPOL (Duke University)
- CIMT-CIP


C. Quality Management System

- GLP Certificate
- ISO 9001:2015



MONITORING VACCINE-SPECIFIC IMMUNE RESPONSES

I. ACUTE INFECTIONS / PROPHYLACTIC VACCINATION




IFN-γ ex-vivo ELISpot

- O/N Resting
- 24h stimulation

Volunteers before/after SARS-CoV-2 or Flu vaccination

COVID-19 convalescent patients from, not vaccinated

II. CHRONIC INFECTIONS / THERAPEUTIC VACCINATION



IFN-γ ELISpot with expansion

- 5 days in vitro expansion

31 HPV+ patients before (pre-vacc) & after vaccination (TP1-3)

2 vaccine-derived peptide pools


IFN-γ ELISpot after 10-days in vitro expansion

- O/N Resting
- 10 days in vitro expansion

2 patients HBV+ before (pre-vacc.) & after (TP1-3) vaccination

5 vaccine-derived peptide pools

III. Solid tumors / personalized neoantigen-based vaccine



IFN-γ ELISpot with 10-days in vitro expansion

- O/N Resting
- 10 days expansion
- 20 peptide pools from neoAg

Patient with solid tumor before/after personalized neoantigen vaccine

RESPONSE CRITERIA – A NEVER ENDING DEBATE

How to define **positivity threshold** (Moodie Z. et al, *Cancer Immunol Immunother.* 2010, 59: 1489):

- Empirical rules or statistical tests (need analysis in quadruplicate)
- Our recommendation** = empirical rule based on both
 - Specific / mock spots ratio > 3-4
 - Specific minus mock spots > 50 spots per million PBMC

