

INTRODUCTION

Biomarker assay development in tissues can pose multiple challenges including complete extraction of the protein from the tissue, wide endogenous ranges, immunocapture (IP) reagent availability, co-interactions with other proteins limiting successful IP extraction. BAG3 is a clinically relevant cardiac biomarker, with mutations linked to familial dilated cardiomyopathy, and exists in complex with its binding partner HSPB8 protein. Simultaneous quantitation of both proteins from limited heart tissue is analytically challenging. Current work showcases the challenges in the development of a multiplexed hybrid LC-MS/MS assay to quantify levels of both BAG3 and HSPB8 from a single human heart tissue sample using complementary offline and online immunoaffinity enrichment strategies. Method development included optimization of tissue homogenization process, identifying the right surrogate matrix, and evaluating multiple enrichment strategies. Upon method optimization, a combined offline and online approach was adopted from a single post-digest sample split, minimizing sample consumption while maintaining assay performance.

CONSTRAINTS AND CONSIDERATIONS


Key Considerations for LC-MS/MS Protein Biomarker Assay Development

When developing an assay for the quantification of protein biomarkers using LC-MS/MS, several factors play a key role in ensuring the assay is robust, reproducible, and fit for its intended context of use. Below is a checklist to consider during **method development**.

1. MATRIX	2. ANALYTE	3. ASSAY CONSIDERATIONS
<ul style="list-style-type: none"> Availability of Matrix: Rare vs readily available Sourcing Healthy vs Diseased Matrix: Consider biological relevance and consistency of source If Tissue Matrix: Homogenization buffer and ratio for consistent recovery 	<ul style="list-style-type: none"> Endogenous: Identification of right surrogate matrix Abundance: Low abundant or high abundant? Need for parallelism and identifying the right diluent Availability of Protein Standard: Recombinant vs cell lysate vs peptide standard Matrix Effect: "Normal levels" vs Disease state levels 	<ul style="list-style-type: none"> Availability of capture reagents: Antibody or affinity reagent availability, specificity, and performance Multiplexing: Basal levels of each biomarker and range identification Establishing Fit-for-Purpose Quality Controls: Consideration to see if 1 set can work for both biomarkers Need for unique peptide to monitor vs any peptide from the biomarker that ionizes well Other assay parameters: Selectivity, Carryover, Dilution integrity, Stability

Careful consideration of these factors during assay development will help ensure a **robust, reproducible, and fit-for-purpose** LC-MS/MS assay for protein biomarker quantification.


HYBRID LC-MS/MS WORKFLOW WITH PEPTIDE-LEVEL CLEAN-UP



METHOD OPTIMIZATION DATA



Parallelism

Parallelism was successfully established up to 200-fold for both biomarkers



Carryover

Single sample was digested and split assay was developed to tackle carryover issue with biomarker 1. A plate-based anti-peptide IP eliminated the need for on-line Ab column and trap before MS analysis. Cause for carryover was identified to be with Ab column and trap for ELL peptide



RESULTS

BIOMARKER 1 PRECISION AND ACCURACY									
Actual Concentration (ng/mL)	Num. Values	Mean	%CV	% Accuracy	Value #1	Value #2	CALIBRATION CURVE		
5	2 of 2	5.14	4.1	102.8	4.99	5.29			
10	2 of 2	9.69	0.3	96.9	9.67	9.71			
20	2 of 2	19.3	3.6	96.7	18.8	19.8			
40	2 of 2	38.4	4.5	96.0	39.6	37.2			
80	2 of 2	81.2	4.3	101.5	78.8	83.7			
160	2 of 2	160	0.0	99.8	160	160			
320	2 of 2	320	1.4	99.9	316	323			
640	2 of 2	663	1.0	103.6	658	668			
1280	2 of 2	1326	3.7	103.6	1360	1291			
2560	2 of 2	2609	4.1	101.9	2684	2534			
4900	2 of 2	4352	3.7	96.7	4467	4238			
5120	2 of 2	5152	0.2	100.6	5146	5157			

Surrogate Curve and Matrix Quality Controls

1x RIPA buffer was used as a surrogate matrix for a 2-plex assay. Based on the endogenous values detected, standard curve range was set from 5-5,120 ng/mL. Due to limited tissue availability, single sample volume of 25 µL was utilized to quantitate levels of both biomarkers implementing a split assay format. A 50x dilution was employed to prepare matrix QCs to bring the quantitation levels within the range of the curve. This decision was supported with the parallelism data obtained.

BIOMARKER 2 PRECISION AND ACCURACY									
Actual Concentration (ng/mL)	Num. Values	Mean	%CV	% Accuracy	Value #1	Value #2	CALIBRATION CURVE		
5	2 of 2	5.17	2.1	103.4	5.09	5.26			
10	2 of 2	9.52	8.2	95.2	10.1	8.97			
20	2 of 2	19.3	3.6	96.3	18.8	19.8			
40	2 of 2	39.0	3.3	97.5	39.9	38.1			
80	2 of 2	80.9	3.0	101.1	79.2	82.6			
160	2 of 2	165	4.9	103.3	159	171			
320	2 of 2	318	3.7	99.4	310	326			
640	2 of 2	667	5.5	104.2	641	693			
1280	2 of 2	1289	1.9	100.7	1271	1307			
2560	2 of 2	2556	1.4	99.9	2581	2531			
4900	2 of 2	4367	2.2	97.1	4435	4300			
5120	2 of 2	5217	3.0	101.9	5105	5329			

CONCLUSION

A multiplexed hybrid LC-MS/MS method was successfully developed and fully validated for a fit for purpose context of use (CoU) format for the simultaneous quantitation of BAG3 and HSPB8 in human heart tissue over a concentration range of 5-5,120 ng/mL with an acceptance criteria set at ±25%. During method development, parameters evaluated included optimal homogenization ratio and buffer selection, surrogate matrix, quality control concentrations, parallelism, selectivity, accuracy and precision, sensitivity, matrix effects, stability, dilution integrity, and reference standard performance. The validated multiplex assay was further qualified in non-human primate and mouse heart and liver tissues, demonstrating translatability across species and matrices. A single-sample, split-workflow multiplex biomarker quantitative hybrid LC-MS/MS method in tissue samples using online and offline IP in tandem was successfully developed and validated under ICH M10 guidance.

Human Heart Tissue		
Biomarker 1		Biomarker 2
460 - 5000 ng/mL	"Normal" Level Ranges	280 - 3000 ng/mL
5-5120 ng/mL	Std Curve Range	5-5120 ng/mL
50x diluted	Quality control basal sample	50x diluted
1x RIPA Buffer	Diluent	1xRIPA Buffer
non-dilute to 200x dilute	Parallelism	non-dilute to 200x dilute
Direct digest with anti-peptide Ab IP	Extraction	Direct digest with anti-peptide Ab IP
Offline anti-peptide Ab IP to minimize Carryover issues	Approach Utilized	On-line Ab column to maximize Recovery

Online Ab column, Trap and Elute Configuration

